




ora
Dental Implant

IMPLANT LINES

INTERNAL HEX



MINI IMPLANT



CONTENT INDEX


BWS® SURFACE	1	OVERVIEW PROSTHETIC	
IMPLANT LINES	2	COMPONENTS - Ø 3.25	19
PACKAGING	3		
<hr/>			
PARALLEL IMPLANT		OVERVIEW PROSTHETIC	
Technical features	4	COMPONENTS - Ø 3.75	22
Reference codes	5		
<hr/>			
ACTIVE IMPLANT		RATCHET	
Technical features	6	CLEANING AND MAINTENANCE	26
Reference codes	7		
<hr/>			
SHORT IMPLANT		MINI IMPLANT	
Technical features	8	Technical features	27
Reference codes	9	Reference codes	28
<hr/>			
SURGICAL TRAY - TRAY "IS"	10	Insertion procedure	29
SURGICAL TRAY - TRAY "IM"	12		
<hr/>			
READING DEPTH NOTCHES AND SHARP DRILLS	14	PRELIMINARY INDICATIONS FOR SURGICAL INSTRUMENT USE	30
DRILL STOP	14		
IMPLANT CONNECTIONS FEATURES AND MEASURES	15	INSTRUMENTS FOR SURGERY WARNING and LEGENDS	31
IMPLANTS INSERTION PROCEDURE	16		
IMPLANT CONNECTION SCREWDRIVERS O-RING REPLACEMENT TOOL	17		
SURGICAL INSTRUMENTS	18		

BWS®


a surface with over 20 years of history

CONSTANT OVER TIME

The capacity of **BWS®** to **retain fibrin**, lets osteoblasts migrate from the bone to the implant surface and reproduce there, **generating new bone** in direct contact with the titanium (contact Osseointegration).



Bone tissue grown in direct contact with the surface **BWS®**




20 µm

SEM HV: 20.00 kV WD: 10.6470 mm
SEM MAG: 4.82 kx Det: SE Detector
View field: 62.05 µm

VEGA\\TESCAN


The process of sandblasting and acid etching the implant surface makes it possible to obtain **optimal values of roughness** creating the strongest fibrin adhesion to the surface and facilitating the bone healing process by **significantly reducing the time.**



2µm

EHT=18.00 kV WD=13 mm Mag=6.50 K X
Photo No.=6159 Detector= SE1

After the surface treatment and the classic washings, Dental Tech implants are additionally cleaned with **Argon Cold Plasma** to minimize carbon contamination. Subsequently, minute controls are performed on the fixture with **scanning electron microscopes (SEM).**



BWS®

- ✓ Packaging in controlled environments
- ✓ Clean room ISO 7
- ✓ Packaging impermeable to micro-organisms
- ✓ Gamma ray sterilisation process guarantee the creation of products that are extremely safe for users and their patients



PARALLEL IMPLANT

Fixture with cylindrical body and a conical apex. Modulating the surgical procedure it is indicated in all bone types; even in the case of non-compact bone it is able to achieve a good primary stability.

You can use it for any type of prosthetic restoration, screwed and cemented. Using the concept of platform switching allows you to better manage the soft tissue in the area of the implant – abutment interface, and reduce peri-implant bone resorption over time.

SHORT IMPLANT

Even if they are 6.0 mm length, allows the realization of surgical procedures without bone graft even in cases of advanced tissue resorption.



ACTIVE IMPLANT

Tapered implant that, thanks to its special spiral design, facilitates the users in the realization of Ridge Expansion procedures. The exceptional self tapping power of the thread, provides an excellent bone condensing and a high primary stability even in very complex clinical cases. Implogic AT is recommended in cases of post extraction implants and in case of poor quality bone.



PACKAGING

ORA Dental Implant GHBH endosseous implants are supplied in sterile packaging which, if undamaged, guarantees the implant is protected from external agents and, if stored correctly, their sterility.



TECHNICAL FEATURES

PARALLEL IMPLANT

BETTER PENETRATION

Spiral profile with hybrid progress: flat and radiating towards the root, triangular-shaped externally, for greater penetration into incompletely prepared sites.

Micro-grooves to limit bone resorption.

The implant's screwing axis can be adjusted.

APICAL DRILLS

Drills with helicoidal progress to enhance stable penetration.

PARALLEL IMPLANT REFERENCE CODES

INTERNAL HEX

Diameter (Ø) mm
Ø 2,20

Length (L) mm	REF
10	FTCP11/SC
11,5	FTCP111/SC
13	FTCP113/SC
16	FTCP117/SC

Recommended surgical sequence



DRLNC DRP200 DRP280 CTK325

Diameter (Ø) mm
Ø 3.75

Length (L) mm	REF
8	I0,3V0A
10	I0,3V10
11,5	I0,3V11
13	I0,3V13
16	I0,3V17

Recommended surgical sequence



DRLNC DRP200 DRP310 DRP325 CTK375
LOW DENSITY HIGH DENSITY

Diameter (Ø) mm
Ø 4.25

Length (L) mm	REF
8	I0,4P0A
10	I0,4P10
11,5	I0,4P11
13	I0,4P13
16	I0,4P17

Recommended surgical sequence



DRLNC DRP200 DRP310 DRP360 DRP375 CTK425
LOW DENSITY HIGH DENSITY

Diameter (Ø) mm
Ø 4.75

Length (L) mm	REF
8	I0,4V0A
10	I0,4V10
11,5	I0,4V11
13	I0,4V13

Recommended surgical sequence



DRLNC DRP200 DRP280 DRP310 DRP360 DRP410 DRP425 CTK475
LOW DENSITY HIGH DENSITY

Diameter (Ø) mm
Ø 0,0-

Length (L) mm	REF
8	I0,000A
10	I0,0010
11,5	I0,0011
13	I0,0013

Recommended surgical sequence



DRLNC DRP200 DRP280 DRP310 DRP360 DRP410 DRP425 CTK475
LOW DENSITY HIGH DENSITY

TECHNICAL FEATURES

ACTIVE IMPLANT

SPIRAL DESIGN

The unusual spiral design simplifies the procedures of Ridge Expansion.

RISK REDUCTION

Less risk of damaging adjacent teeth and perforation of the lingual and/or buccal cortical plates.

SELF-TAPPING COIL


Exceptional self-tapping capability which provides improved bone condensation and increased primary stability, even in highly complex clinical cases.

BONE MAINTENANCE OVER TIME

Allows a greater reduction of bone osteotomy to be achieved, which results in lower bone loss and reduced surgical trauma.

OPTIMAL CHOICE OF POSITIONING

Allows a change in direction in order to achieve the optimum position of restoration, especially in post-extraction sites.



ACTIVE IMPLANT REFERENCE CODES

INTERNAL HEX




Diameter (Ø) mm
Ø 3.75

Length (L) mm	REF
8	11,3V08
10	11,3V10
11,5	11,3V11
13	11,3V13
16	11,3V16




Diameter (Ø) mm
Ø 4.25

Length (L) mm	REF
8	11,4P08
10	11,4P10
11,5	11,4P11
13	11,4P13
16	11,4P16



Diameter (Ø) mm
Ø 4.75

Length (L) mm	REF
8	11,4V08
10	11,4V10
11,5	11,4V11
13	11,4V13



TECHNICAL FEATURES

SHORT IMPLANT

SELF-TAPPING COIL

Self-tapping coil with double principle thread for increased contact with the bone and greater primary stability.

BONE MAINTENANCE OVER TIME

Polished coronal chamfer and implant collar are designed to better manage the biological width and maintain the level of bone over time.

SPIRE GEOMETRY

The geometry of the spire aids osseous healing, both qualitatively and quantitatively.


IMPROVED PENETRATION

Four wide cutting zones for greater penetration capacity and to gather bone fragments, therefore reducing compression.

FACILITATES POSITIONING THE DEVICE IN THE SURGICAL SITE


Tapered apical portion to facilitate centring of the device in the surgical site, even in cases of under preparation due to poor bone density, or to achieve greater primary stability.

ATRAUMATIC APEX



SHORT IMPLANT REFERENCE CODES

INTERNAL HEX



Diameter (Ø) mm
Ø 4.25

Length (L) mm	REF
7	I .,EP.7

Recommended surgical sequence

Drill	DRLNC	DRP200	DRP280	DRP310	DRP360
Stop		STC2506	STC2506	STC2506	STC3406

Diameter (Ø) mm
Ø 4.75

Length (L) mm	REF
7	I .,EV.7

Recommended surgical sequence

Drill	DRLNC	DRP200	DRP280	DRP310	DRP360	DRP410
Stop		STC2506	STC2506	STC2506	STC3406	STC3406

Diameter (Ø) mm
Ø 0,0•

Length (L) mm	REF
7	I .,00.7

Recommended surgical sequence


Drill	DRLNC	DRP200	DRP280	DRP310	DRP360	DRP410	DRP425
Stop		STC2506	STC2506	STC2506	STC3406	STC3406	STC3406

SURGICAL TRAY - "TRAY IS"


REF TRAY IS

DIMENSIONS

142x104 mm - h 61 mm



SURGICAL TRAY - "TRAY IS"



Ø3.75
L 6
P...E




L 18,5
TW...IL

L 12,5
TW...IC



(at choice) (at choice)



COMPARTMENT INSIDE
DRILL STOP

L 14,3
GCG...P

L 8,3
GCG...E



(at choice) (at choice)

L 9
K10A9




ERS

SPARE

SCREW
DRIVER C.

D.EXTENSION



Ø3.6 Ø3.75 Ø4.1 Ø4.25

Ø4.25 Ø4.75

Ø 4.25
L 9,5
CTKEP0



Ø 4.75
L 9,5
CTKEV0



Ø 3.6
L 23
DRP11



Ø 3.75
L 23
DRP10



Ø 4.1
L 23
DRP11



Ø 4.25
L 23
DRPE0




SURGICAL TRAY - "TRAY IM"


REF TRAY IM

DIMENSIONS


176x143 mm - h 63,5 mm



SURGICAL TRAY - "TRAY IM"




P...I



COMPARTMENT INSIDE
DRILL STOP



L 9
K10A9





ERS



SPARE

SCREW
DRIVER C.

D.EXTENSION



READING DEPTH NOTCHES AND SHARP DRILLS




COUNTERSINK



DRILL STOP

Diameter (Ø) mm **Stop**
Ø ϵ ,0



Length (L) mm	REF
6	STCP0-1
7	STCP0-V
8	STCP0-Λ
10	STCP01-
11,5	STCP01 I

Length (L) mm **Parallel Drill**
L $\pi\pi$

Diameter (Ø) mm	REF
2.0	DRPπ-.
2.8	DRPπΛ.
3.1	DRPπI.
3.25	DRPππ0

Diameter (Ø) mm **Stop**
Ø 0,0



Length (L) mm	REF
6	STCπε-1
7	STCπε-V
8	STCπε-Λ
10	STCπε1-
11,5	STCπε1 I

Length (L) mm **Parallel Drill**
L $\pi\pi$

Diameter (Ø) mm	REF
3.6	DRPπ1.
3.75	DRPπV0
4.1	DRPε1-
4.25	DRPεπ0

Countersink

Diameter (Ø) mm	REF
3.25	CTKππ0

Countersink

Diameter (Ø) mm	REF
4.25	CTKεπ0

Countersink

Diameter (Ø) mm	REF
3.75	CTKπV0



Countersink

Diameter (Ø) mm	REF
4.75	CTKεV0




IMPLANT CONNECTIONS


FEATURES AND MEASURES



IMPLANT HANDWEEL CONNECTION




IMPLANT CA CONNECTION



IMPLANT RATCHET CONNECTION


	IMPLANT HANDWEEL CONNECTION	IMPLANT CA CONNECTION	IMPLANT RATCHET CONNECTION
Ø 3,20 - L 11,0	REF CM1370	REF IC0370	REF IM0370
Ø 3,20 - L 7,0 SHORT	REF 3,20,07	REF CAS-ZIM	REF 3,20,08
Ø 3,20 - L 11,0 LONG	REF 3,20,11	REF CAL-ZIM	REF 3,20,10

Ø Diameter mm - L Length mm




IMPLANT HANDWEEL CONNECTION

Allows removal of the implant from the ampoule and the start of insertion in the surgical site.



IMPLANT CA CONNECTION

Allows removal of the implant from the ampoule and its insertion in the surgical site using the contra-angle screwdriver.



IMPLANT RATCHET CONNECTION

A tool to be connected to the ratchet to complete insertion of the implant. It does not permit removal as it does not have an O-Ring seal.

IMPLANTS INSERTION PROCEDURE

WITH MANUAL CONTRA-ANGLE IMPLANT CONNECTION

Insert the direct manual contra-angle screwdriver into the implant with a slight rotating motion to allow the correct coupling of the two hexagons (implant - screwdriver) and remove the implant. (Fig. 5)

Begin insertion of the implant in the alveolar surgery (Fig. 6) after having set the following parameters on the surgical unit:

- 1) Bi-phase procedure (submerged) RPM 20-15
Torque max. 40-35 Ncm
- 2) Monophasic procedure realized with submerged implants and healing screws, with deferred load RPM 20-15
Torque max. 45-40 Ncm
- 3) Monophasic procedure with immediate load/prosthesis RPM 20-15
Torque is incremental from 20 to 70 Ncm

If a surgical unit with good torque control is available, both in quantity and quality, it is possible to terminate insertion of the implant with the contra-angle; if the opposite is true, insert the device in the alveolar surgery as long as the power of the machine permits and complete the insertion manually proceeding as follows:

IMPLANT RATCHET CONNECTION

Ensure that the tool is inserted in the position suitable for screwing and turn until the implant reaches the desired position. (Fig. 7)

Complete the insertion of the implant using the dynamometric wrench connected to the direct screwdriver of the REF. IMA20.005 / 20.004 / 325 ratchets. At times it is necessary to use the extensions, short REF. PMC115 and long REF. 110026 to connect to the tools described above. (Fig. 8)



Fig.5



Fig.6



Fig.7



Fig.8

IMPLANT CONNECTION SCREWDRIVERS

O-RING REPLACEMENT TOOL

DESCRIPTION

The tools (DRO 375) facilitates replacement of the O-Ring on the screwdrivers.



Warning The O-Ring support tool is made of PMMA and, therefore, it cannot be sterilised in an autoclave.




Fig.11









Fig.12



Fig.13



RATCHET
REF **1110I**

<p>BALL ABUTMENTS ADAPTOR</p> <p>for ball abutment Ø 2.25</p>  <p>inox</p> <table border="1"> <tr><td>REF</td></tr> <tr><td>RDS110</td></tr> </table>	REF	RDS110	<p>SCREWDRIVERS ADAPTOR</p>  <p>inox</p> <table border="1"> <tr><td>REF</td></tr> <tr><td>Short TW0001C</td></tr> <tr><td>Long TW0001L</td></tr> </table>	REF	Short TW0001C	Long TW0001L	<p>OVERDENTURE ABUTMENT ADAPTOR</p>  <p>Ti₅</p> <table border="1"> <tr><td>REF</td></tr> <tr><td>ADL150</td></tr> </table>	REF	ADL150								
REF																	
RDS110																	
REF																	
Short TW0001C																	
Long TW0001L																	
REF																	
ADL150																	
<p>MUA ADAPTOR</p>  <p>inox</p> <table border="1"> <tr><td>REF</td></tr> <tr><td>TW001A</td></tr> </table>	REF	TW001A	<p>ISO ADAPTOR</p>  <p>inox</p> <table border="1"> <tr><td>Length (L) mm</td><td>REF</td></tr> <tr><td>7</td><td>IS011V</td></tr> </table>	Length (L) mm	REF	7	IS011V	<p>HEX SCREWDRIVER</p>  <p>inox</p> <table border="1"> <tr><td>Length (L) mm</td><td></td><td>REF</td></tr> <tr><td>8,8</td><td>Short</td><td>GCG0024</td></tr> <tr><td>14,8</td><td>Long</td><td>GCG0030</td></tr> </table>	Length (L) mm		REF	8,8	Short	GCG0024	14,8	Long	GCG0030
REF																	
TW001A																	
Length (L) mm	REF																
7	IS011V																
Length (L) mm		REF															
8,8	Short	GCG0024															
14,8	Long	GCG0030															

SURGICAL INSTRUMENTS

DYNAMOMETRIC RATCHET



inox

Length (L) mm	REF
	••1101

ADAPTOR FOR DYNAMOMETRIC RATCHET



inox

ISO connection for ratchet

Length (L) mm	REF
7	ISO•V•

EXTENSION



inox

Length (L) mm	REF
12,5	11••27

EXTENSION FOR DRILL




inox

Length (L) mm	REF
9	K10A9


SCREWDRIVERS

inox




Length (L) mm	REF
4,5	GMX1••
11,5	GMM•0•
18	••110P

HEX CA DRIVER



inox

Length (L) mm	REF
8,3	Short GCG••2E
14,3	Long GCG••2•




DEPTH GAUGE

Length (L) mm	REF
	••11E•

OVERVIEW PROSTHETIC COMPONENTS

DIAMETER 3.25



H Height mm
 Ø Diameter mm
 ML Laser Marking

COMPONENTS FOR IMPRESSIONS AND MODELS

*Fastening screw included

PICK UP IMPRESSION COPY



Ø	REF
3.8	TPR ² ...

TRANSFER IMPRESSION COPING



REF
TST ² 0

IMPLANT ANALOGUES



REF
AGL ² 1 ²

IMPLANTS



REF
FTC ² 11/SC
FTC ² 11/SC
FTC ² 12/SC
FTC ² 17/SC

HEALING ABUTMENT



H	REF
2	VG ² 02
4	VG ² 0E
6	VG ² 01

COMPONENTS FOR CEMENTED PROSTHETICS



REF VFD²2²
 Fastening screw included and available as a replacement

STRAIGHT ABUTMENT

Ti₀




H	Ø	ML	REF
1,5	3.8	325	MAS ² 10
3	3.8	325	MAS ² 20

Ø	REF
3.8	PLT ² 0

20Ncm Torque adapter REF TW0001

ANGLED ABUTMENT

Ti₀



H	Ø	ML	REF
1,5	3.8	15/325°	MPG ² 11
1,5	3.8	25/325°	MPG ² 21
3	3.8	15/325°	MPG ² 12
3	3.8	25/325°	MPG ² 22

20Ncm Torque adapter REF TW0001


CASTABLE ABUTMENT

Pmma



REF
PCA ² 0

COMPONENTS FOR PROSTHETICS SCREWED AT THE IMPLANT LEVEL



REF VFD²2²
 Fastening screw included and available as a replacement

TEMPORARY CYLINDER

PEEK



ML	REF
325E	PKR ² ...
325R	PKR ² ...

20Ncm Torque adapter REF TW0001

CYLINDER ABUTMENT

Ti₀



ML	REF
325E	PPE ² 0
325R	PPR ² ...

20Ncm Torque adapter REF TW0001

CASTABLE ABUTMENT

Pmma



REF
PCA ² 0
PCR ² ...

20Ncm Torque adapter REF TW0001

OVERVIEW PROSTHETIC COMPONENTS

DIAMETER 3.25

COMPONENTS FOR PROSTHETICS SCREWED TO AN ABUTMENT



REF **VBTRP..**

* Included and available as a replacement

ABUTMENT MUA

Ti₀

H



H	REF
1	BTRPI..
2,5	BTRPPO

20Ncm Torque adapter
REF TW0001

CYLINDER

Ti₀



***VBTRP..**

REF
CITRP..

Pmma



REF
CBTRP..

20Ncm Torque adapter
REF TW0001

IMPRESSION TRANSFER

Ti₀



***VBTRP..**

REF
TBTRP..

ABUTMENT ANALOGUE

Ti₀



REF
ABTRP..

PROTECTIVE CAP



Ti₀

***VBTRP..**

REF
GBTRP..

CAD-CAM COMPONENTS

* Fastening screw included and available as a replacement

SCAN ABUTMENT

Ti₀



***VFXRPPO**

REF
SCANRPPO

Digital CAD-CAM Intraoral Scan and Laboratory Scan.
For single cemented and screwed elements.
For multiple cemented elements.

DIGITAL ANALOGUE

Ti₀



REF
AGLRPDG

Analog for digital models, specific for applications through the manufacture of models made with 3D printing/prototyping. The characteristic shape with rounded edges, allows easy insertion into the model seat, without interference and friction with the resinous material of the models.
The apical screw allows to always obtain a total working stability. This prosthetic component must be used through the Libraries.

TI BASE SIRONA

Ti₀



***VFXRPPO**

H	REF	REF
0,5	● PSSRPPO	● PSSRPOR
1	● PSSRPI..	● PSSRPI..R
2	● PSSRPP..	● PSSRPP..R

Digital CAD-CAM bonding technique.
For single cemented and screwed elements.
For multiple cemented elements.

20Ncm Torque adapter
REF TW0001

MUA BONDING BASE

Ti₀



***VBTRP..**

REF
BCMRPO

Digital CAD-CAM and traditional bonding technique.
For multiple elements screwed into an MUA Pillar.

20Ncm Torque adapter
REF TW0001

CAD-CAM COMPONENTS

SCAN BTA

Ti₀



REF

SCANBTA

Screw included

Suitable for digital CAD-CAM technique for intraoral and laboratory scans. For multiple screw-retained elements.

DIGITAL ANALOGUE

Ti₀



REF

ABT··DG**

Analogue for digital models, specific for applications through the manufacture of models made with 3D printing/prototyping. The characteristic shape with rounded edges allows an easy insertion in the seat model made, without interference and friction with the resinous material of the models. The apical screw allows to always obtain a total working stability. This prosthetic component must be used in conjunction with the Libraries.

OVERDENTURE COMPONENTS - BALL ATTACHMENT

BALL ABUTMENT O-RING

Pack. 10 pcs

Ti₅



REF

POR225

REF

ORG225

BALL ABUTMENT ANALOG

Ti₅



REF

AAF0**

BALL ABUTMENT

Ti₅




H	REF
0,5	ASF**··
1,5	ASF**10
3	ASF**·

20Ncm Torque adapter REF RDS225

OVERVIEW PROSTHETIC COMPONENTS

DIAMETER 3.75



H Height mm
Ø Diameter mm
ML Laser Marking

IMPLANTS



REF	REF	REF	REF
10,3V0A	10,4P0A	10,4V0A	10,000A
10,3V10	10,4P10	10,4V10	10,0010
10,3V11	10,4P11	10,4V11	10,0011
10,3V13	10,4P13	10,4V13	10,0013
10,3V17	10,4P17		



REF	REF	REF
11,3V0A	11,4P0A	11,4V0A
11,3V10	11,4P10	11,4V10
11,3V11	11,4P11	11,4V11
11,3V13	11,4P13	11,4V13
11,3V17	11,4P17	



REF	REF	REF
10,4P07	10,4V07	10,0007

HEALING ABUTMENT



Parallel



Taper



Taper Large

H	REF	H	REF	H	REF
2	00,01	2	00,04	2	00,0V
4	00,02	4	00,00	4	00,0A
6	00,03	6	00,07	6	00,09

COMPONENTS FOR IMPRESSIONS AND MODELS

* Fastening screw included

PICK UP IMPRESSION COPY

Ti ₀	Ø	ML	REF
	4.5	Z	00,11
	5.6	ZL	00,12
Ti ₀	Ø	ML	REF
	3.5	Z	00,10

*00,0A

TRANSFER IMPRESSION COPING

Ti ₀	ML	REF
	Z	00,14
	Z	00,13
Ti ₀	ML	REF
	Z	00,13

*00,0A1

IMPLANT ANALOGUES

Ti ₀	ML	REF
	Z	00,10

COMPONENTS FOR CEMENTED PROSTHETICS



REF 00,0A2
 Fastening screw included and available as a replacement

STRAIGHT ABUTMENT

Ti ₀	Ø	ML	REF
	3.75	Z	00,17

H	Ø	ML	REF
1,5	4.5	Z	00,1V
3	4.5	Z	00,1A
1,5	5.5	ZL	00,19
3	5.5	ZL	00,20

2 Ncm Torque adapter REF TW0001

CASTABLE ABUTMENT



Pmma

ML	REF
ZE	00,29

2 Ncm Torque adapter REF TW0001


OVERVIEW PROSTHETIC COMPONENTS

DIAMETER 3.75

COMPONENTS FOR CEMENTED PROSTHETICS

ANGLED ABUTMENT

Ti₀



H	Ø	ML	REF
1,5	4.5	Z15	•0,•21
1,5	4.5	Z25	•0,•23
3	4.5	Z15	•0,•22
3	4.5	Z25	•0,•24
1,5	5.5	ZL15	•0,•20
1,5	5.5	ZL25	•0,•2V
3	5.5	ZL15	•0,•27
3	5.5	ZL25	•0,•28

20Ncm Torque adapter REF TW0001



REF •0,•A2
Fastening screw included and available as a replacement

COMPONENTS FOR PROSTHETICS SCREWED AT THE IMPLANT LEVEL

TEMPORARY CYLINDER

Peek



ML	REF
ZE	•0,•22
ZR	•0,•21

20Ncm Torque adapter REF TW0001

CYLINDER ABUTMENT

Ti₀



ML	REF
ZE	•0,•23
ZR	•0,•24

20Ncm Torque adapter REF TW0001



REF •0,•A2
Fastening screw included and available as a replacement

CASTABLE ABUTMENT

Pmma




ML	REF
ZE	•0,•29
ZR	•0,•20

20Ncm Torque adapter REF TW0001

COMPONENTS FOR PROSTHETICS SCREWED TO AN ABUTMENT

MUA STRAIGHT ABUTMENT

Ti₀



H	Ø	REF
1	4.5	•0,•20
3	4.5	•0,•21

Supplied with transfer handle

20Ncm Torque adapter REF TW0080

MUA ANGLED ABUTMENT

Ti₀



*REF •0,•A2

H	Ø	ML	REF
1	4.5	Z17	•0,•2V
3	4.5	Z17	•0,•28
1	4.5	Z30	•0,•29
3	4.5	Z30	•0,•20

Supplied with transport and parallelization screw*

20Ncm Torque adapter REF TW0001

MUA PROTECTION CAP

Peek



ML	REF
ZM	•0,•21

Use only on MUA abutments

Pack. 2 pcs

MUA CYLINDER

Ti₀



ML	REF
ZM	•0,•22

Use only on MUA abutments

Lock manually

Pmma



ML	REF
ZM	•0,•23

MUA PROSTHETIC SCREW

Ti₀



REF
•0,•24

Use only on MUA abutments

Lock manually

Ti₀



REF
•0,•20

Pack. 2 pcs

OVERVIEW PROSTHETIC COMPONENTS

DIAMETER 3.75

COMPONENTS FOR PROSTHETICS SCREWED TO AN ABUTMENT

* Fastening screw included

MUA PRECISION TRANSFER



Ti₀ *•0,•ΛΕ

ML	REF
ZM	•0,•ΕΤ

MUA ABUTMENT ANALOGUE



Ti₀

ML	REF
ZM	•0,•ΕV

MUA SCAN ABUTMENT



Ti₀

ML	REF
ZM	•0,•0V

Screw included

Suitable for digital CAD-CAM technique for intraoral and laboratory scans. For multiple screw-retained elements.

MUA DIGITAL ANALOGUE



Ti₀

ML	REF
ZD	•0,•0Λ

Analogue for digital models, specific for applications through the manufacture of models made with 3D printing/prototyping. The characteristic shape with rounded edges allows an easy insertion in the seat model made, without interference and friction with the resinous material of the models. The apical screw allows to always obtain a total working stability. This prosthetic component must be used in conjunction with the Libraries.

CAD-CAM COMPONENTS

* Fastening screw included and available as a replacement

SCAN ABUTMENT




Ti₀

*•0,•Λ0

ML	REF
Z	•0,•00

Digital CAD-CAM Intraoral Scan and Laboratory Scan. For single cemented and screwed elements. For multiple cemented elements.

DIGITAL ANALOGUE



Ti₀

ML	REF
ZD	•0,•0Γ

Analogue for digital models, specific for applications through the manufacture of models made with 3D printing/prototyping. The characteristic shape with rounded edges, allows easy insertion into the model seat, without interference and friction with the resinous material of the models. The apical screw allows to always obtain a total working stability. This prosthetic component must be used through the Libraries.

TI BASE SIRONA



H *•0,•Λ0

Ti₀

H	ML	REF	REF
0,5	Z	•0,•Εϱ	•0,•0Γ
1	Z	•0,•0•	•0,•0Γ
2	Z	•0,•01	•0,•0Ε

Digital CAD-CAM and traditional bonding technique. For single cemented and screwed elements. For multiple cemented elements.

20Ncm Torque adapter
REF TW0001

MUA BONDING BASE



Ti₀

*•0,•ΕΕ

ML	REF
ZM	•0,•ΕΛ


Digital CAD-CAM and traditional bonding technique. For multiple screwed elements on MUA pillar.

20Ncm Torque adapter
REF TW0001

OVERVIEW PROSTHETIC COMPONENTS


DIAMETER 3.75

DIGITAL ANALOGUE - INDICATIONS OF USE



OVERDENTURE COMPONENTS – BALL ATTACHMENT

BALL ABUTMENT



H	REF
0,5	•0,•09
1,5	•0,•1•
3	•0,•11
5	•0,•11


20Ncm Torque adapter REF RDS225

BALL ABUTMENT O-RING



Pack. 10 pcs


BALL ABUTMENT ANALOG



RATCHET CLEANING AND MAINTENANCE



RATCHET
REF 20,001



The dynamometric ratchet, after each use, must be disassembled for cleaning. This maintenance operation does not require any tools.


Completely unscrew the screw **(A)**, remove the whole pawl **(B)** and then the flexible dynamometric bar **(C)**. Once disassembled, clean according to the instructions for use and maintenance attached to the device, brush with non-metallic rigid bristles, even in hollow areas with pipe cleaner for a complete removal of biological residues.

Once the cleaning and disinfection phase has been completed, reassemble the ratchet using the reverse disassembly procedure, making sure to match the pin **(D)** in the housing dedicated.

MINI IMPLANT

The MINI IMPLANT system meets the growing clinical need to have small diameter implants for instant stabilisation of total prostheses. Designed for long-term rehabilitation and conceived for excellent clinical results.

EXCEPTIONALLY EASY
Implant characteristics make the surgical phase very easy. The ergonomics of supplied components facilitate prosthetic procedures. Hence, implants can be inserted and the prosthesis can be stabilised in just one session.



EXCELLENT RESISTANCE

The implant is a monocomponent made of Titanium Gr5 for maximum mechanical resistance.

SMALL PROFILE

The diameter (barely 2.7 mm) allows to place the implant in the thin crestal bone to avoid bone regeneration procedures.

MAXIMUM BONE SURFACE CONTACT


The development of implant macrotopography and the clinically tested surface obtained with the BWS[®] system ensures excellent primary stability of the device and a high BIC (Bone Implant Contact).

MINIMALLY INVASIVE SURGERY

The dentist can choose whether to insert the implant with a traditional or flapless technique.

REFERENCE CODES

MINI IMPLANT




Sphere diameter (Ø) mm
Ø P, P0

Lenght (L) mm	REF
10	SPHPVI ·/S
11,5	SPHPVI I/S
13	SPHPVI P/S

Ti₀

Load direction



Static load

Breakage at N 1500

Stress resistance

N 505 x 5,000,000 cycles

No breakages


The tests were performed on devices with Ø2mm to evaluate either the same or higher diameters.



INITIAL CYLINDRICAL DRILL

inox

Lenght (L) mm	REF
2,0	DRPP ··



TAP

Ti₀

Lenght (L) mm	REF
L 10 - LI 5	MSPHPV ··



HANDWHEEL

Ti₀

Lenght (L) mm	REF
6	AMC · I1



IMPLANTS SCREWDRIVERS

inox

Lenght (L) mm	REF
10	RDSrro

RATCHET

inox

REF
P ·, ·· I



ABUTMENT

Ti₀



REF
PTSPH

BALL ABUTMENT O-RING

Ti₅



REF
POR225



REF
ORG225

Pack. 10 pcs

BALL ABUTMENT ANALOG

Ti₅



REF
AAF225

Remove the device Mini Implant, which is connected to the plastic cap, from the ampoule by concurrently pulling and gently rotating the cap. (Fig. 14)



Fig.14

Carry the implant into the mouth with the cap/support, and use it to start placing the implant in the osteotomy site. Screw the implant onto the bone until it reaches the stability level that allows to extract the support from the device by pulling upwards. (Fig. 15)



Fig.15

Complete the insertion of the Mini Implant by using manual key assembled with the dedicated adaptor to screw it on, leaving the entire hexagonal portion that is under the sphere outside the soft tissue. This will prevent the O-RING retention device from causing compression of soft tissues. (Fig. 16)



Fig.16

Other instruments can also be used as an alternative to the pawl. (Fig. 17)




Fig.17

Insertion of the implant must ensure that O-RING retainers are correctly in place. Hence the need to ensure a distance of at least 7mm between the osteotomies. (Fig. 18)




Fig.18



Warning Do not exceed 20 rpm and 55 Ncm of torque when screwing on the implant.

PRELIMINARY INDICATIONS FOR SURGICAL INSTRUMENT USE

PREVENTION

Besides correct and continuous long-term maintenance, wear and tear of the instruments can also be prevented and slowed down.

In the first place every instrument must only be used for the envisaged and indicated use.

The instruments used must be cleaned immediately after the end of surgery.

Remove residue and encrustations only with soft brushes and NOT with metal brushes.

When envisaged, disassemble the instruments and deeply clean the cavity. The devices must be fully immersed in the most appropriate detergents or disinfectants for the material, and left to rest for a period of time that never exceeds the manufacturer's instructions.

After disinfecting them, rinse thoroughly with water and dry the devices with a clean and dry cloth. Complete with a jet of compressed air.

CLEANING PHASES

PACKAGING AND STERILITY

- ORA Implant tools are supplied as non sterile in heat-sealed Pouches in containing the leaflet.
- ORA Implant tools can be used again and therefore it has to be washed and sterilised prior to their usage.

Dental Tech validated the following cleansing and disinfection method:

MANUAL CLEANING

- Just after the use of ORA Implant equipment, place the equipment into a container with a peracetic acid based solution at concentration of 2% (NO GLUTARALDEHYDE OR SODIUM HYPOCHLORITE), as long as 18 minutes.
- After-ward rinse carefully.

MANUAL DISINFECTION

- Place the equipment into a container with a peracetic acid based solution at concentration of 4% (NO GLUTARALDEHYDE OR SODIUM HYPOCHLORITE), as long as 15 minutes.
- Rinse generously
- Examine the equipment and make sure there are no organic remains. Carefully scrub the outer parts with a non-metal bristled brush.

MANUAL RINSE

- Place the equipment into ultrasound bath, and wash it for approx. 18 minute and then rinse carefully.

DRY

- Perfectly dry the equipment, seal it individually with material suitable for moist heat sterilisation.

CHECK

After the cleaning phases, check that none of the instruments presents signs of corrosion, contamination or damage. Especially use a magnifying lens to check the most concealed areas, the joints and the handles.

If any contamination is detected, repeat the cleaning procedure.

In case of damage, dispose of the instrument as established by the laws in force for waste management.

STERILISATION

Sterilise in a steam autoclave saturated with distilled water by using a systematically validated and controlled sterilisation method, according to provisions laid down by standard ISO 1:2007-17665 "Sterilisation of healthcare products" (as amended). Requirements for validation and routine control of moist heat sterilisation in healthcare facilities".

- Dental Tech validated the following Autoclave moist heat sterilization cycle:

3 minutes

134 °C



Warning The use of suitable protection during cleaning and sterilisation of contaminated instruments enhances personal safety during these phases.

Since Dental Tech tools are manufactured in different materials, they shall be washed and sterilized one by one.

PRESERVATION

After the sterilisation phase, the instruments must be preserved in the sterilised package in a dry, dust-free place, far from heat sources. The bags must only be opened before use.

The storage period of sterilised items must not exceed the period recommended and indicated on the bag.

DISPOSAL PROCEDURES

At the end of its life the medical device must be disposed of according to the methods established by national laws in force for waste management.

INSTRUMENTS FOR SURGERY WARNINGS AND LEGENDS

INSTRUMENT FOR SURGERY

The surgical instrumentation of the Dental Tech Implant System is simple and essential, responding to every clinical need and treatment protocol. All drills and components are laser marked, to allow preparation of the implant site correctly to the established depth, and a predictable and safe positioning of the implant. The instruments are available individually or in sets with different types of surgical kit.

HOW TO USE THE SURGICAL INSTRUMENTS

So as not to cause mechanical and/or thermal damage to bone tissue in the zone in which the implant is to be inserted, and to obtain a congruous surgical site (indispensable to achieving good osseointegration of the implant) some fundamental rules must be respected:

- Use drills with gradual diameter progression: the same instruments must not be used for more than 25 osteotomies;
- Do not exceed 800 RPM during the osteotomy;
- Do not exceed 20 RPM in the event of tapping with the contra-angle;
- Ensure, during the osteotomy, that the instruments work in axis;
- Do not exert lateral pressure during the osteotomy and tapping;
- The osteotomy must be performed exercising light pressure and back and forth movements on the axis of the instrument;
- Use generous irrigation with physiological solution, both during drilling and tapping of the surgical site;
- Ensure that during the intervention the irrigation canals of the instruments are clear;
- Avoid categorically, during surgery, the cooling of instruments and the implant site with the air-water syringes tips.

NON-ROTATING INSTRUMENT

The non-rotating instrument is compatible with all Dental Tech implant systems.







WARNINGS

RESPONSABILITY The use of non-original components, produced by third-parties may compromise the functionality of the implants and their elements, compromising the final result and voiding the guarantee of the manufacturer. The application of the product occurs outside the control of Dental Tech and is the sole responsibility of the end user. We accept no liability for any damage resulting from such activities.











INSTRUCTIONS FOR USE These are to be considered solely as recommendations. This information is not sufficient and does not exempt the user from ensuring the adequacy of the product for its intended use through continued training.

VALIDITY This nullifies all previous versions. The images, the content and the products illustrated are subject to modification without warning.

MATERIALS LEGEND

	Gold Alloy
	Surgical Stainless Steel
	Polyetereeterechetone
	Polymethylmethacrylate
	Titanium gr.V ELI for medical use
	Polymer

PACKAGING SYMBOLS LEGEND

	Lot number
	Sterilized by gamma rays
	Not sterile
	Product code
	Reusable
	Use by
	Non-reusable
	Attention, consult the supplied documentation
	Directive 94/93/CEE conformity mark
	Notified body identification

BWS® are registered trademarks of DentalTech * S.r.l.
Any reproduction or publication is permitted only with the written authorization of DentalTech * S.r.l.

IMPLANT LINE

IMPLASSIC FTP



CONTENT INDEX

BWS® SURFACE	1	INSTRUMENTS	18
TECHNICAL FEATURES AND PACKAGING	2	PRELIMINARY INDICATIONS FOR SURGICAL INSTRUMENT USE	19
SURGICAL PROCEDURE AND REFERENCE CODES	3	INSTRUMENTS FOR SURGERY WARNING and LEGENDS	20
READING DEPTH NOTCHES and SHARP DRILLS	4	BIBLIOGRAPHY	21
STOP INSERTION AND REMOVAL PROCEDURE	5		
SURGICAL TRAY	6		
SCREWDRIVER Features and Measures	8		
IMPLANTS INSERTION PROCEDURE	9		
SURGICAL INSTRUMENTS	10		
HEALING ABUTMENT PROSTHETIC CONNECTION	11		
COMPONENTS FOR PROSTHETIC CONOMETRY	12		
OVERVIEW PROSTHETIC COMPONENTS	13		
PROSTHETIC DIGITAL COMPONENTS	16		
OVERVIEW PROSTHETIC COMPONENTS FOR TORONTO BRIDGE, SCREWED BRIDGE AND FULL ARCH	17		




BWS®



a surface with over 20 years of history

CONSTANT OVER TIME

The capacity of **BWS®** to **retain fibrin**, lets osteoblasts migrate from the bone to the implant surface and reproduce there, **generating new bone** in direct contact with the titanium (contact Osseointegration).



Bone tissue grown in direct contact with the surface **BWS®**



20 µm

SEM HV: 20.00 kV WD: 10.6470 mm
SEM MAG: 4.82 kx Det: SE Detector
View field: 62.05 µm

VEGA\\TESCAN Dental Tech

The process of sandblasting and acid etching the implant surface makes it possible to obtain **optimal values of roughness** creating the strongest fibrin adhesion to the surface and facilitating the bone healing process by **significantly reducing the time**.



2µm

EHT=18.00 kV WD=13 mm Mag=6.50 K X
Photo No.=6159 Detector= SE1

After the surface treatment and the classic washings, Dental Tech implants are additionally cleaned with **Argon Cold Plasma** to minimize carbon contamination. Subsequently, minute controls are performed on the fixture with **scanning electron microscopes (SEM)**.



BWS®

- ✓ Packaging in controlled environments
- ✓ Clean room ISO 7
- ✓ Packaging impermeable to micro-organisms
- ✓ Gamma ray sterilisation process guarantee the creation of products that are extremely safe for users and their patients


TECHNICAL FEATURES



Conometric connection at 6°, with hexagonal position index and screw through, extremely precise and stable.



Smooth collar 0,75mm. The eccentric course between implant and connection diameter offers an anatomical path to the prosthetic component.



The geometric peculiarity of the cortical spiral allows to obtain an high primary stability, even in the presence of a few millimeters of bone.



Thanks to the flat shape of the central loop, the FTP implant allows the condensation of the bone matrix during the insertion of the fixture.



Apical spiral with progressive course allow greater directionality in insertion, in addition to the high primary stability in poor quality bone.




The atraumatic apex, without cutting areas, makes the implant suitable even in cases where it is necessary to safeguard anatomical structures, such as maxillary sinus and alveolar nerve.

PACKAGING

ORA Dental Implant GHBH endosseous implants are supplied in sterile packaging which, if undamaged, guarantees the implant is protected from external agents and, if stored correctly, their sterility.



Protective implant cap supported by a titanium ring. (Surgical colour code)



Transparent ampoule

Protective closure screw cap. (Prosthetic colour code)

SURGICAL PROCEDURE AND REFERENCE CODES



Diameter (Ø) mm $\text{Ø } r,vo$



Lenght (L) mm	REF
10	FTPr,vo /SC
11,5	FTPr,vo I /SC
13	FTPr,vo II /SC
16	FTPr,vo III /SC

Prosthetic colour code ●

Surgical colour code ●

Cover screw included

Recommended surgical sequence



Diameter (Ø) mm $\text{Ø } \epsilon,ro$



Lenght (L) mm	REF
8	FTPϵ,ro /SC
10	FTPϵ,ro I /SC
11,5	FTPϵ,ro II /SC
13	FTPϵ,ro III /SC
16	FTPϵ,ro IV /SC

Prosthetic colour code ●

Surgical colour code ●

Cover screw included

Recommended surgical sequence



Diameter (Ø) mm $\text{Ø } \epsilon,vo$


Lenght (L) mm	REF
8	FTPϵ,vo /SC
10	FTPϵ,vo I /SC
11,5	FTPϵ,vo II /SC
13	FTPϵ,vo III /SC

Prosthetic colour code ●

Surgical colour code ●

Cover screw included

Recommended surgical sequence




★ It is recommended if the cortical bone is very persistent.



Warning All DRP drills are 0,8 mm longer than the implant. In the planning stage and while drilling in proximity to vital anatomical structures, this added length must be considered.

READING DEPTH NOTCHES AND SHARP DRILLS



DRILL STOP

Ø 0,0
L mm implant

Diameter (Ø) mm Ø 0,0	Stop
Lenght (L) mm	REF
6	STC ϵ 0·1
7	STC ϵ 0·V
8	STC ϵ 0·A
10	STC ϵ 0I·
11,5	STC ϵ 0II
13	STC ϵ 0I ϵ
16	STC ϵ 0I7

3.1
L

Lenght (L) mm L ϵ ϵ	Parallel Drill
Diameter (Ø) mm	REF
2.0	DRP ϵ ··
2.3	DRP ϵ · ϵ
2.8	DRP ϵ A·
3.1	DRP ϵ I·
3.25	DRP ϵ ·0

Ø 0,0
L mm implant

Diameter (Ø) mm Ø 0,0	Stop
Lenght (L) mm	REF
6	STC ϵ ·1
7	STC ϵ ·V
8	STC ϵ ·A
10	STC ϵ I·
11,5	STC ϵ II
13	STC ϵ I ϵ
16	STC ϵ I7

4.1
L

Lenght (L) mm L ϵ ϵ	Parallel Drill
Diameter (Ø) mm	REF
3.6	DRP ϵ ·1·
3.75	DRP ϵ V0
4.1	DRP ϵ I·
4.25	DRP ϵ ·0

STOP INSERTION AND REMOVAL PROCEDURE

STOP INSERTION

Hold the drill by the stalk and insert the stop, with the retentive flaps facing towards the drill, until it comes into contact with the metal stop located on the drill itself. (Fig. 3 - 2 - 1)

STOP REMOVAL

Hold the stop and remove the drill, pulling on the side of the stalk.



CORRECT



Fig.1



Fig.2



Fig.3



WRONG

STOP WRONG INSERTION

The Stop insertion with the wings facing the tip of the drill is incorrect. (Fig. 4)



Fig.4

DEPTH STOP FOR DIFFERENT LENGTHS

ADVANTAGES


- Optimum control of depth during Preparation of the surgical site, even in conditions of poor visibility in the operative field;
- Reduction of surgical risk;
- Reduction of operator stress;
- Greater patient safety;
- Facilitates the insertion and removal of the drill stop and increased safety during surgery for the doctor and assistant, the cutting portion of the instrument is never touched by the operators.

SURGICAL TRAY - "TRAY IM"


REF TRAY IM

DIMENSIONS


176x143 mm - h 63,5 mm



COMPARTMENT INSIDE
DYNAMOMETRIC RATCHET P, . . . I



SURGICAL TRAY - "TRAY M"



(at choice)



(at choice)



(at choice)



(at choice) (at choice)



(at choice) (at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice) (at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)



(at choice)




(at choice)




(at choice)

SCREWDRIVERS

FEATURES AND MEASURES



- L 7,0
- L 11
- L 10,0



- L 8
- L 12
- L 17,0

L Length mm



DIRECT MANUAL CONTRAANGLE SCREWDRIVER

Hexagonal head for ratchet



FINAL STEP SCREWDRIVER

It allows the implant removal from the ampoule and its insertion into the surgical site by ratchet or handwheel.



SCREWDRIVERS'
References in millimeters.
Useful for the management of hard and soft tissues in the positioning of the implants.

WITH MANUAL SCREWDRIVER

Insert the screwdriver (IMAFTP001-IMAFTP002-IMAFTP003), connected to the handwheel (AMC016), into the implant making a slight rotation to allow good matching of the two hexagons (implant - screwdriver) and remove the implant. (Fig. 1)

Begin insertion of the implant in the alveolar surgical site using the manual screwdriver. Where bone density permits, it is possible complete insertion of the implant using the manual wrenches. (Fig. 2)



Fig.1



Fig.2

To remove, exercise a slight lateral movement, right and left, in order to free the conometric matching.

WITH DIRECT CONTRA-ANGLE SCREWDRIVER

Insert the direct manual contra-angle screwdriver into the implant with a slight rotating motion to allow the correct coupling of the two hexagons (implant - screwdriver) and remove the implant. (Fig. 3)

Begin insertion of the implant in the alveolar surgery (Fig. 4) after having set the following parameters on the surgical unit:

- 1) Bi-phase procedure (submerged) RPM 20-15 Torque max. 40-35 Ncm
- 2) Monophasic procedure realized with submerged implants and healing screws, with deferred load RPM 20-15 Torque max. 45-40 Ncm
- 3) Monophasic procedure with immediate load/prosthesis RPM 20-15 Torque is incremental from 20 to 70 Ncm

If a surgical unit with good torque control is available, both in quantity and quality, it is possible to terminate insertion of the implant with the contra-angle; if the opposite is true, insert the device in the alveolar surgery as long as the power of the machine permits and complete the insertion manually proceeding as follows.



Fig.3



Fig.4

To remove, exercise a slight lateral movement, right and left, in order to free the conometric matching.

FINAL SCREWDRIVER

Ensure that the tool is inserted in the position suitable for screwing and turn until the implant reaches the desired position. (Fig. 5)

Complete the insertion of the implant using the dynamometric wrench connected to the direct screwdriver of the ratchets. At times it is necessary to use the extensions, short REF. PMC115 and long REF. 110026 to connect to the tools described above. (Fig. 6)




Fig.5



Fig.6


To remove, exercise a slight lateral movement, right and left, in order to free the conometric matching.




PROSTHETIC EXTRACTOR

Lenght (L) mm	REF
21	EXPFTP
28	EXPFTPL

HANDWHEEL



Lenght (L) mm	REF
6	AMC·17



DYNAMOMETRIC RATCHET

REF
CCD·V·

ADAPTOR FOR DYNAMOMETRIC RATCHET



ISO connection for ratchet

Lenght (L) mm	REF
7	ISORV·



PARALLEL PIN

REF
CPTV·EV

RATCHET

REF
PGII··




EXTENSION

Lenght (L) mm	REF
7	PMCI10
12,5	II··P7




EXTENSION FOR DRILL

Lenght (L) mm	REF
9	KI0A9




HEX SCREWDRIVER

Lenght (L) mm		REF
4,5	Micro	GMXI··
11,5	Extra Short	GMMPO·
13,5	Long	··II0P



HEX CA DRIVER

Lenght (L) mm		REF
8,3	Short	GCG··PE
14,3	Long	GCG··P·



DEPTH GAUGE

REF
··IIE·

HEALING ABUTMENT PROSTHETIC CONNECTION



ORA Dental Implant GHBH's FTP implant line offers clinicians **versatility of use** that makes this type of implant suitable for any surgical indication.

The **1° conometric connection**, with hexagonal position index and through screw, allows an accurate and stable matching of the prosthetic components.

Conometric matching at 6° between fixture and abutment, with the presence of a hexagonal index to facilitate the positioning of the abutment.



IMPLASSIC FTP

Ø mm	Length mm
3.75	16 - 13 - 11,5 - 10
4.25	16 - 13 - 11,5 - 10 - 8
4.75	13 - 11,5 - 10 - 8




Important Warning

Excessive torques can compromise the hexagonal shape of the screws and screwing tools, causing impediments, even irreversible, during operating and prosthetic phases. The recommended tightening torques for the screws are summarized in the following table:

SCREW DESCRIPTION	INSTRUMENT	TORQUE Ncm
Surgical Screw	Manual screwdrivers	manually 10/8Ncm
Healing Abutment	Manual screwdrivers	manually 10/8Ncm
Transfer Screw	Manual screwdrivers	manually 10/8Ncm
Fixing Screw Abutment MUA (M1,4)	Manual screwdrivers	manually 10/8Ncm
Scan Abutment screws	Manual screwdrivers	manually 10/8Ncm
Fixing Screw Abutment	Adaptor for dynamometric ratchet Contra-Angle Screwdriver	20Ncm



Given the importance of tightening torque, it is recommended to always monitor the perfect functionality of the tightening tools, evaluating carefully the tools and subjecting them to constant maintenance. It is always recommended to start tightening the screws using manual screwdrivers and, only for the determination of the correct tightening torque, for screws that have a specific torque, use the appropriate tools to impress the indicated torque.



ANATOMIC HEALING ABUTMENT

Height (HT) mm	H1	H2	REF
4	2	2	VGFTPε•0•
6	3	3	VGFTP1•0•

Ti₀



CYLINDRICAL HEALING ABUTMENT

Height (H) mm	REF
4	VGFTPρ0ε•
6	VGFTPρ01•

Ti₀

COMPONENTS FOR PROSTHETIC CONOMETRY

FASTENING SCREW

Included and available as a replacement

REF **VMTP·P·V**
Tighten to 20 Ncm

Ti₀



CONOMETRIC ABUTMENT 0°



(HT) mm	REF
1	MTPPVI·
2	MTPPV·

CONOMETRIC ABUTMENT 10°



(HT) mm	REF
1	PTPPVI I
1	PTPPVI IR
2	PTPPVI P
2	PTPPVI PR

CONOMETRIC ABUTMENT 20°



(HT) mm	REF
1	PTPPV I
1	PTPPV IR
2	PTPPV P
2	PTPPV PR

CONOMETRIC CAP FOR WELDING



Ti₀

REF
CPDε··

CONOMETRIC CAP FOR AESTHETIC MATERIAL



Ti₀

REF
CPP··ε

CONOMETRIC TRANSFER SNAP



Pmma

REF
CAI·00

CONOMETRIC ANALOG



Ti₀

REF
AFTPCAP

HEALING ABUTMENTS PROSTHETIC CONOMETRY ABUTMENTS



Pmma

REF
CGDI··

OVERVIEW PROSTHETIC COMPONENTS

FASTENING SCREW
Included and available
as a replacement



REF **VMTP·P_V**
Tighten to 20 Ncm

Ti₀

STRAIGHT ABUTMENT SNAP



Ti₀

(HT) mm	REF
2	PDFTP1··
3	PDFTP2··

ANGLED ABUTMENT 1° SNAP



Ti₀

(HT) mm	REF
2	PAFTP1··1
3	PAFTP1··2

ANGLED ABUTMENT 2° SNAP



Ti₀



(HT) mm	REF
2	PAFTP2··1
3	PAFTP2··2

THREE-PART PRECISION TRANSFER FOR PICK-UP TECHNIQUE



Ti₀

REF
TPRFTP



Transfer for Pick Up technique. Used with a perforated impression tray, it allows the removal of the central body of the Transfer by extracting the anti-rotation hexagonal connection, in order to facilitate the removal of the impression, in the event of disparallelisms between implants.

Indications Analog alignment:

the flat side (A) is always corresponding to the connection hexagon side (B). For a practical alignment it is recommended to keep the flat side (A) and connection hexagon side (B) in correspondence with the external body flat side (C).

PLASTER ANALOG



Ti₀

REF
AGFTP_PV

TRANSFER SNAP

For PDFTP and
PAFTP abutments



Peek

REF
SNAPFTP

Pack. 4 pcs

OVERVIEW PROSTHETIC COMPONENTS

EQUATOR FTP



(HT) mm	REF
0,0	·P·FTPE·0
1	·P·FTPE1
2	·P·FTPE2
ε	·P·FTPEε

2 · Ncm Torque adapter

KIT SMART BOX



REF
330SBC

CAPS ASSORTMENT KIT



REF
19PECE

VIOLET CAP GR 2700 (strong)



REF
1ε·CEV

Pack. 4 pcs

CLEAR CAP GR 1800 (standard)



REF
1ε·CET

Pack. 4 pcs

PINK CAP GR 1200 (soft)



REF
1ε·CER

Pack. 4 pcs

YELLOW CAP GR 600 (extra-soft)



REF
1ε·CEG

Pack. 4 pcs

BLACK CAP FOR LABORATORY



REF
1ε·CEN

Pack. 4 pcs

STAINLESS STEEL HOUSING



REF
1ε1CAE

Pack. 4 pcs

SMARTBOX HOUSING WITH BLACK CAP



REF
33·SBE

LABORATORY ANALOG



REF
1εεAE

Pack. 2 pcs

IMPRESSION COPINGS NORMAL SIZE – OT EQUATOR



REF
·εεCAIN

Pack. 2 pcs

TRANSPARENT PROTECTIVE DISKS



REF
1··PD

Pack. 10 pcs

OVERVIEW PROSTHETIC COMPONENTS

FASTENING SCREW

Included and available as a replacement

REF **VMTP·P**
Tighten to 20 Ncm



Ti₀

ABUTMENT TYPE SIRONA S



HT

Ti₀

(HT) mm	REF	
0,0	PSFTP·0·	⬡
I	PSFTP I · ·	⬡
P	PSFTP P · ·	⬡

Suitable for digital and traditional bonding technique CAD-CAM. For single screwed elements on the implant. **Portion coronal compatible SIRONA.**

ROTATING ABUTMENT TYPE SIRONA S



HT

Ti₀

(HT) mm	REF	
0,0	PSFTP·0·R	●
I	PSFTP I · ·R	●
P	PSFTP P · ·R	●

Suitable for digital and traditional bonding technique CAD-CAM. For single screwed elements on the implant. **Portion coronal compatible SIRONA.**

SCAN ABUTMENT SIRONA



Plastica

REF
7EPIPII

Pack. 36 pcs

Digital CAD-CAM Intraoral Scan and Laboratory Scan. For single cemented and screwed elements. On SIRONA pillar.

CYLINDER



Ti₀

REF	
CTFTP&P·	⬡
CTFTP&P·R	●

CASTABLE ABUTMENT



Pmma

REF	
CCFTP&V·	⬡
CCFTP&V·R	●

PREMILLED



*

Ti₀

REF
PRMLI I FTP
PRMLI 7 FTP

* FASTENING SCREW

Included and available as a replacement

REF **VMTP·P**
Tighten to 20 Ncm

The Pre-Milled are indicated to design and carry out customized abutment processing depending on the patient's clinical case. The Dental Technician, through the ORA Libraries, has the possibility to determine the morphology of the abutment with the main CAD Designers. The Pre-Milled are available in 2 diameters (11.5mm and 16mm) for each prosthetic platform, depending on the inclination designed and are compatible with the most widespread attachment-machine on the market.

PROSTHETIC DIGITAL COMPONENTS

SCAN ABUTMENT



* **SCAN ABUTMENT SCREW** included

Also available as a replacement.
REF **VFSFTP**

Ti₀

REF
SCANFTP

Digital CAD-CAM Intraoral Scan and Laboratory Scan.
For single cemented and screwed elements - multiple cemented elements.

DIGITAL ANALOG



Ti₀

REF
AGFTPVDG

Analog for digital models, specific for applications through the manufacture of models made with 3D printing/prototyping. The characteristic shape with rounded edges, allows easy insertion into the model seat, without interference and friction with the resinous material of the models. The apical screw allows to always obtain a total working stability. This prosthetic component must be used through the ORA Libraries.

DIGITAL ANALOG - INDICATIONS OF USE



CORRECT POSITIONING



WARNING
DO NOT orient the Scan Abutment in the unsuitable and aligned secondary position

It is necessary to match up the smaller portion of the Scan Abutment, which is always oriented on the hexagonal side of the connection, with the side of the external square shape of the analogous digital body.

OVERVIEW PROSTHETIC COMPONENTS FOR TORONTO BRIDGE, SCREWED BRIDGE AND FULL ARCH

ABUTMENT MUA TRANSFER



Ti₀

REF
HTMEI · 7

ABUTMENT MUA ANALOG



Ti₀

REF
HLM · · EI

PROTECTION CAP MUA



Peek

REF
HPMEI · ·

Pack. 2 pcs
Prosthetic screw NOT included

STRAIGHT ABUTMENT MUA



HT

Ti₀

(HT) mm	REF
1,0	PEFTP10 ·
2	PEFTP2 · ·

ANGLED ABUTMENT 1V° MUA



HT

Ti₀

(HT) mm	REF
2	PEFTP1V2 ·
3	PEFTP1V3 ·

ANGLED ABUTMENT 2° MUA



HT

Ti₀

(HT) mm	REF
2	PEFTP2 · 2 ·
3	PEFTP2 · 3 ·

CYLINDER MUA



Ti₀

REF
HMT · · EI

Prosthetic screw NOT included

CASTABLE ABUTMENT MUA



Pmma

REF
HMC EI · ·

Prosthetic screw NOT included

PROSTHETIC SCREW



Ti₀

REF	
VPCEM	Short
VPLEM	Long

Pack. 2 pcs

ABUTMENT MUA DIGITAL COMPONENTS

SCAN MUA



MA

Ti₀

REF
SCANMA

Screw included (REF VPCEM)

Suitable for digital CAD-CAM technique, for intraoral and laboratory scans. For multiple screwed elements.

DIGITAL ANALOG MUA



MA

Ti₀

REF
HLM · · EI DG

Analog for digital models, specific for applications through the manufacture of models made with 3D printing/prototyping. The characteristic shape with rounded edges, allows easy insertion into the model seat, without interference and friction with the resinous material of the models. The apical screw allows to always obtain a total working stability. This prosthetic component must be used through the ORA Libraries.

BONDING BASE FOR ABUTMENT MUA



Ti₀

REF
BCMHEX


Suitable for digital CAD-CAM technique, for intraoral and laboratory scans. For multiple screwed elements.

INSTRUMENTS



DYNAMOMETRIC RATCHET
REF **CCD·V·**

CLEANING AND MAINTENANCE



The dynamometric ratchet, after each use, must be disassembled for cleaning. This maintenance operation does not require any tools.

Completely unscrew the screw **(A)**, remove the whole pawl **(B)** and then the flexible dynamometric bar **(C)**. Once disassembled, clean according to the instructions for use and maintenance attached to the device, brush with non-metallic rigid bristles, even in hollow areas with pipe cleaner for a complete removal of biological residues.

Once the cleaning and disinfection phase has been completed, reassemble the ratchet using the reverse disassembly procedure, making sure to match the pin **(D)** in the housing dedicated.

TIGHTENING TOOLS FOR DYNAMOMETRIC RATCHET

ADAPTOR FOR FASTENING SCREWS

inox



	REF
Short	TW0001C
Long	TW0001L

OCTA ABUTMENT ADAPTOR

inox



REF
TW··Λ·

ISO CONNECTION FOR RATCHET


inox



Lenght (L) mm	REF
7	ISOPV·

HEX SCREWDRIVER

inox



Lenght (L) mm		REF
8,8	Short	GCG0024
14,8	Long	GCG0030

PRELIMINARY INDICATIONS FOR SURGICAL INSTRUMENT USE

PREVENTION

Besides correct and continuous long-term maintenance, wear and tear of the instruments can also be prevented and slowed down.

In the first place every instrument must only be used for the envisaged and indicated use.

The instruments used must be cleaned immediately after the end of surgery.

Remove residue and encrustations only with soft brushes and NOT with metal brushes.

When envisaged, disassemble the instruments and deeply clean the cavity. The devices must be fully immersed in the most appropriate detergents or disinfectants for the material, and left to rest for a period of time that never exceeds the manufacturer's instructions.

After disinfecting them, rinse thoroughly with water and dry the devices with a clean and dry cloth. Complete with a jet of compressed air.

PACKAGING AND STERILITY

- ORA Implant tools are supplied as non sterile in heat-sealed Pouches in containing the leaflet.
- ORA Implant tools can be used again and therefore it has to be washed and sterilised prior to their usage.

Dental Tech validated the following cleansing and disinfection method:

MANUAL CLEANING

- Just after the use of ORA Implant equipment, place the equipment into a container with a peracetic acid based solution at concentration of 2% (NO GLUTARALDEHYDE OR SODIUM HYPOCHLORITE), as long as 18 minutes.
- After-ward rinse carefully.

MANUAL DISINFECTION

- Place the equipment into a container with a peracetic acid based solution at concentration of 4% (NO GLUTARALDEHYDE OR SODIUM HYPOCHLORITE), as long as 15 minutes.
- Rinse generously
- Examine the equipment and make sure there are no organic remains. Carefully scrub the outer parts with a non-metal bristled brush.

MANUAL RINSE

- Place the equipment into ultrasound bath, and wash it for approx. 18 minute and then rinse carefully.

DRY

- Perfectly dry the equipment, seal it individually with material suitable for moist heat sterilisation .

CHECK

After the cleaning phases, check that none of the instruments presents signs of corrosion, contamination or damage. Especially use a magnifying lens to check the most concealed areas, the joints and the handles.

If any contamination is detected, repeat the cleaning procedure.

In case of damage, dispose of the instrument as established by the laws in force for waste management.

STERILISATION

Sterilise in a steam autoclave saturated with distilled water by using a systematically validated and controlled sterilisation method, according to provisions laid down by standard ISO 1:2007-17665 "Sterilisation of healthcare products" (as amended). Requirements for validation and routine control of moist heat sterilisation in healthcare facilities".

- Dental Tech validated the following Autoclave moist heat sterilization cycle:

3 minutes

134 °C



Warning The use of suitable protection during cleaning and sterilisation of contaminated instruments enhances personal safety during these phases.

Since ORA tools are manufactured in different materials, they shall be washed and sterilized one by one.

PRESERVATION

After the sterilisation phase, the instruments must be preserved in the sterilised package in a dry, dust-free place, far from heat sources. The bags must only be opened before use.

The storage period of sterilised items must not exceed the period recommended and indicated on the bag.

DISPOSAL PROCEDURES

At the end of its life the medical device must be disposed of according to the methods established by national laws in force for waste management.

INSTRUMENTS FOR SURGERY WARNINGS AND LEGENDS

INSTRUMENT FOR SURGERY

The surgical instrumentation of the Dental Tech Implant System is simple and essential, responding to every clinical need and treatment protocol. All drills and components are laser marked, to allow preparation of the implant site correctly to the established depth, and a predictable and safe positioning of the implant. The instruments are available individually or in sets with different types of surgical kit.

HOW TO USE THE SURGICAL INSTRUMENTS

So as not to cause mechanical and/or thermal damage to bone tissue in the zone in which the implant is to be inserted, and to obtain a congruous surgical site (indispensable to achieving good osseointegration of the implant) some fundamental rules must be respected:

- Use drills with gradual diameter progression: the same instruments must not be used for more than 25 osteotomies;
- Do not exceed 800 RPM during the osteotomy;
- Do not exceed 20 RPM in the event of tapping with the contra-angle;
- Ensure, during the osteotomy, that the instruments work in axis;
- Do not exert lateral pressure during the osteotomy and tapping;
- The osteotomy must be performed exercising light pressure and back and forth movements on the axis of the instrument;
- Use generous irrigation with physiological solution, both during drilling and tapping of the surgical site;
- Ensure that during the intervention the irrigation canals of the instruments are clear;
- Avoid categorically, during surgery, the cooling of instruments and the implant site with the air-water syringes tips.

NON-ROTATING INSTRUMENT

The non-rotating instrument is compatible with all ORA implant systems.







WARNINGS

RESPONSABILITY The use of non-original components, produced by third-parties may compromise the functionality of the implants and their elements, compromising the final result and voiding the guarantee of the manufacturer. The application of the product occurs outside the control of ORA and is the sole responsibility of the end user. We accept no liability for any damage resulting from such activities.











INSTRUCTIONS FOR USE These are to be considered solely as recommendations. This information is not sufficient and does not exempt the user from ensuring the adequacy of the product for its intended use through continued training.

VALIDITY This nullifies all previous versions. The images, the content and the products illustrated are subject to modification without warning.

MATERIALS LEGEND

 Au	Gold Alloy
 inox	Surgical Stainless Steel
 Peek	Polyetereeterechetone
 Pmma	Polymethylmethacrylate
 Ti ₆	Titanium gr.V ELI for medical use
 Plastic	Polymer

PACKAGING SYMBOLS LEGEND

 LOT	Lot number
 STERILE R	Sterilized by gamma rays
 NON STERILE	Not sterile
 REF	Product code
 RIUTILIZZABILE	Reusable
	Use by
	Non-reusable
	Attention, consult the supplied documentation
	Directive 94/93/CEE conformity mark
	Notified body identification
0123	

BWS® are registered trademarks of DentalTech * S.r.l. Any reproduction or publication is permitted only with the written authorization of DentalTech * S.r.l.

BIBLIOGRAPHY

Abrahamsson I, Zitzmann NU, Berglundh T, Wennerberg A, Lindhe J. Bone and soft tissue integration to titanium implants with different surface topography: an experimental study in the dog. *Int.J.Oral Maxillofac Implants* 2001; 32-323:(3)16.

Abrahamsson I, Zitzmann NU, Berglundh T, Linder E., Wennerberg A, Lindhe J. The mucosal attachment to titanium implants with different surface characteristics: an experimental study in dogs. *J Clin Periodontal* (5)29 ;2002: 55-448.

The Role of Surface Topography
Herman, *J Perio* 1130-68:1117;1997.

Micro-threads eliminate bone loss due to crestal disuse atrophy
Hansson, *Clin Oral Imp Res*, 1999.

Topografia della neoformazione ossea perimplantare: studio sperimentale
G. Petrone, G. Iezzi, M. Piattelli, A. Scarano
Dipartimento di scienze Odontostomatologiche,
Università "G. D'Annunzio" Chieti- Pescara.

Surface Chemistry Effects of topographic Modification of Titanium Dental Implant Surfaces:
1. Surface Analysis
M. Morra, dr. chem / C. Cassinelli, dr. Biol / G. Bruzzone, MD
A. Capri, MD / G. Di Santi, MD / R. Giardino, MD / M. Fini, MD.
Int. JOMI 45-18:40 ;2003

Surface Chemistry Effects of topographic Modification of Titanium Dental Implant Surfaces:
2. In Vitro Experiments
M. Morra, dr. chem / C. Cassinelli, dr. Biol / G. Bruzzone, MD
A. Capri, MD / G. Di Santi, MD / R. Giardino, MD / M. Fini, MD.
Int. JOMI 52-18:46 ;2003

Valutazione della precisione della connessione tra moncone ed impianto
Benedicenti S.* / Balboni C.** / Maspero F. * / Benedicenti A.*
Quintessence International 4/3 bis 2001

Adesione cellulare epiteliale su superfici di titanio sabbiate e acidificate: studio in vitro
I. Vozza / A. Scarano* / S. Rossi / M. Quaranta
Supplemento n.1 a *Doctor OS* anno XIV n.1 gennaio 2003

Valutazione istologica della risposta ossea a una nuova superficie implantare sabbiate e mordenzate: uno studio sperimentale sul coniglio
Antonio Scarano / Giovanna Iezzi* / Alessandro Quaranta** / Adriano Piattelli*
Implantologia orale numero 2 marzo 2007

Dentista moderno ottobre 2011
Progettazione e realizzazione di una superficie implantare dalla decontaminazione all'osteointegrazione
Chiara Giamberini / Angelo Tagliabue / Dino Azzalin / Giorgio Santarelli

Int J Periodontics Restorative Dent. 2006 Feb; 17-9 :(1)26
Platform switching: a new concept in implant dentistry for controlling postrestorative crestal bone levels.
Lazzara RJ / Porter SS.

I Vela-Nebot X, et al.
Benefits of an implant platform modification technique to reduce crestal bone resorption.
Implant Dent 320-15:313;2006



ORA DENTAL IMPLANT GMBH

Headquarters :
Mühlenstraße 8a
14167 Berlin, Germany

ORA DENTAL IMPLANT GMBH

Headquarters :
Mühlenstraße 8a
14167 Berlin, Germany



DENTAL TECH SRL

Via G. Di Vittorio, 12/10
20826 Misinto (MB), Italy