




ora
Dental Implant

IMPLANT LINES

INTERNAL HEX



MINI IMPLANT

ACTIVE IMPLANT

SPIRAL DESIGN

The unusual spiral design simplifies the procedures of Ridge Expansion.

RISK REDUCTION
Less risk of damaging adjacent teeth and perforation of the lingual and/or buccal cortical plates.


SELF-TAPPING COIL
Exceptional self-tapping capability which provides improved bone condensation and increased primary stability, even in highly complex clinical cases.

BONE MAINTENANCE OVER TIME
Allows a greater reduction of bone osteotomy to be achieved, which results in lower bone loss and reduced surgical trauma.

OPTIMAL CHOICE OF POSITIONING
Allows a change in direction in order to achieve the optimum position of restoration, especially in post-extraction sites.




ACTIVE IMPLANT REFERENCE CODES INTERNAL HEX




Diameter (Ø) mm
Ø 3.75

Length (L) mm	REF
8	11.3708
10	11.3710
11,5	11.3711
13	11.3713
16	11.3716



Diameter (Ø) mm
Ø 4.25

Length (L) mm	REF
8	11.4208
10	11.4210
11,5	11.4211
13	11.4213
16	11.4216



Diameter (Ø) mm
Ø 4.75

Length (L) mm	REF
8	11.4708
10	11.4710
11,5	11.4711
13	11.4713

