



ora
Dental Implant

IMPLANT LINES

INTERNAL HEX



MINI IMPLANT

PARALLEL IMPLANT

BETTER PENETRATION
Spiral profile with hybrid progress: flat and radiating towards the root, triangular-shaped externally, for greater penetration into incompletely prepared sites.

Micro-grooves to limit bone resorption.

The implant's screwing axis can be adjusted.

APICAL DRILLS
Drills with helicoidal progress to enhance stable penetration.

PARALLEL IMPLANT REFERENCE CODES INTERNAL HEX

Diameter (Ø) mm
Ø 3.25

Length (L) mm	REF
10	FTC3210/SC
11,5	FTC3211/SC
13	FTC3213/SC
16	FTC3216/SC

Recommended surgical sequence



Diameter (Ø) mm
Ø 3.75

Length (L) mm	REF
8	10.3708
10	10.3710
11,5	10.3711
13	10.3713
16	10.3716

Recommended surgical sequence



Diameter (Ø) mm
Ø 4.25

Length (L) mm	REF
8	10.4208
10	10.4210
11,5	10.4211
13	10.4213
16	10.4216

Recommended surgical sequence



Diameter (Ø) mm
Ø 4.75

Length (L) mm	REF
8	10.4708
10	10.4710
11,5	10.4711
13	10.4713

Recommended surgical sequence



Diameter (Ø) mm
Ø 5.50

Length (L) mm	REF
8	10.5508
10	10.5510
11,5	10.5511
13	10.5513

Recommended surgical sequence

